

BMS Lead Barriers



Nuclear Medicine Barrier

Specifically designed shorter and wider for shielding nuclear medicine personnel from patient-emitted radiation.

	Dimensions	Lead Equiv.
Window Pane:	36" x 24"	2mm
Opaque Panel:	36" x 34"	2mm

Cat# 56-602 Overall: 36"W x 61"H **Call**



Regular Barrier

This rugged mobile barrier can easily fit into the most confining of spaces and provides exceptional visibility.

	Dimensions	Lead Equiv.
Window Pane:	30" x 24"	2mm
Opaque Panel:	30" x 38"	2mm

Cat# 56-603 Overall: 30"W x 75"H **Call**



Jumbo Barrier

Large 30" x 60" shatter resistant CLEAR-Pb lead plastic window provides superb full height visibility.

	Dimensions	Lead Equiv.
Window Pane:	30" x 60"	2mm
Opaque Panel:	30" x 12"	-

Cat# 56-604 Overall: 30"W x 75"H **Call**



Porta Shield

Mobile radiation protection shield covered with a stain resistant, washable nylon fabric. Available in three lengths (24" shown). Pole adjusts from 36" - 60" in height. 24"W Base.

Choose 0.5mm or 0.35mm Lead Equivalency.

Small Cat# PTS1 24" x 24" Shield **\$470**

Medium Cat# PTS2 24" x 36" Shield **\$551**

Large Cat# PTS3 24" x 48" Shield **\$625**

More **Barriers** on back page!



Burlington Medical Supplies, Inc.

BMS Lead Barriers



Deluxe Barrier

The large 48" wide shatter resistant CLEAR-Pb lead plastic window provides unrestricted visibility and maximum protection for two or more people.

	Dimensions	Lead Equiv.
Window Pane:	48" x 36"	2mm
Opaque Panel:	48" x 36"	2mm

Cat# 56-616 Overall: 48"W x 75"H **Call**

Ultra Barrier

Provides full-body protection to groups of medical personnel assisting in or observing a procedure.

	Dimensions	Lead Equiv.
Window Pane:	72" x 36"	2mm
Opaque Panel:	72" x 36"	2mm

Cat# 56-618 Overall: 72"W x 75"H **Call**

Order Today!

Call 1-800-221-3466

Contact your local BMS Sales Rep
or Email info@burmed.com



Burlington Medical Supplies, Inc.

3 Elmhurst Street • PO Box 3194 • Newport News, VA 23603

©2011 - Information and pricing were correct at time of printing and are subject to change without notice. 11/11