

BMS Laser Protective Eyewear

Two popular filters available in three frame styles!



808 Black



808 Silver



33 Fitover

AKG-5 Homium/Yag/Co2 Filter*

Wavelength	Optical Density	Color	Visible Light Transmission
900-1000nm	5+	Clear	Excellent Visibility
1000-2400nm	7+		
2900-10600nm	7+		

Frame		
808 Black	808 Silver	33 Fitover
LS-KG5-808B	LS-KG5-808S	LS-KG5-33
\$210	\$210	\$200



808 Black



808 Silver



33 Fitover

Multiwave YAG, Alexandrite Diode Filter*

Wavelength	Optical Density	Color	Visible Light Transmission
710-1200nm	7+	Light Blue	64%

Frame		
808 Black	808 Silver	33 Fitover
LS-BG42-808B	LS-BG42-808S	LS-BG42-33
\$250	\$250	\$240

*Customer must determine proper application and use of BMS laser safety products. BMS is in no way liable for any damages arising out of the use, inability to use or incorrect use of our laser safety products.

Visit www.burmed.com

or Call 1-800-221-3466 to order!



Burlington Medical Supplies, Inc.

3 Elmhurst Street • PO Box 3194 • Newport News, VA 23603

©2009 - Information and pricing were correct at time of printing and are subject to change without notice. 07/09